Observation

Microscope

Observation

Microscope

**Ref :
571 042**


Français – p 1

English – p 7

Version : 7003

Microscope polarisant binoculaire

Binocular polarizing microscope

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
NOTES

1 Technical specifications

- Sturdy and stable stand.
- Binocular head **(D)** rotating 360° with observation tubes inclined 30°.
- Dioptic adjustment on one of the observation tubes **(F)**.
- Inter-pupillary distance adjustable from 55 to 75 mm.
- Analyser inbuilt in the head, retractable using tab **(A)**.
- Retractable polarizer rotating 360° **(B)**.
- Eyepieces: wide field 10 x **(E)**.
- 4-position turret fitted with the following objectives **(G)**:

Achromatic objectives	4 x	10 x	40 x
Digital aperture	0.1	0.25	0.65
mounting	fixed	fixed	retractable

- Rapid focusing using coarse and fine-focusing knob **(J)**. Tightness of coarse-focusing knob adjustable using ring
- Circular (Ø 120 mm) rotating stage with 2/10^{ths} graduation over 360° **(H)** with 2 brad-head centring screws and 2 slide holders
- 1.25 digital aperture Abbe condenser with rack and pinion control, iris diaphragm and retractable filter holder **(I)** (blue filter supplied).
- Inbuilt lighting in the base with on/off switch and light intensity adjuster, 6 Volts -20 Watt halogen light.
- Weight of microscope: 5.5 kg.

2 Setting the microscope up


1. Take the microscope out of its box and remove the various protective packaging elements. Keep the packaging for long-term storage or for transport.
2. Fit the binocular head and lock the base using the knurled screw.
3. Screw in the three X 4, X 10 and X 40 objectives onto the turret in increasing order of magnification.
4. Place the oculars in the ocular tubes.
5. The microscope can be used in non-polarized light (e.g. for histology) by removing the polarizing filter and working without the analyser filter; to do this push in the tab **(A)** (on the right side of the head).
6. To use the microscope in polarized light, a very dark background must be obtained (maximum extinction); to do this place the polarizer **(B)** on the light condenser **(C)**, fit the analyser by pulling out the tab **(A)** (on the right side of the head), and rotate the polarizing filter until maximal extinction is obtained. The economy polarizing microscope is equipped with very intense lighting far in excess of that needed for working in polarized light. This intense light source is particularly useful for work in standard light (e.g. in histology, non-polarized light, etc.). To obtain maximal extinction in polarized light (almost black to the eye) all that is required is to reduce the intensity of the light. Maximum extinction is obtained when the polarizing filter and the analysing filter are in crossed positions (90° to each other). The microscope is now ready to be used in polarized light.

3 Using and adjusting the microscope

To move the microscope lift it only by the curved upper part of the stand, placing the other hand under the stand.

- Turn on the on/off switch and rotate the light intensity knob on the base to obtain the required lighting.
- Put the thin slide in place.
- The observation head which can be rotated over 360° and the tubes inclined at a 30° angle provide excellent ergonomics and make prolonged use comfortable.
- The parfocal lenses make it easier to focus when changing magnifications. The X 40 objective is removable in order to avoid damage to the preparation of slide in case of accidental contact.
- Focusing is done by moving the stage vertically using the focusing knob for coarse focusing and a coaxial micrometer knob for fine focusing. The ring on the base of the coarse-focusing knob makes it possible to adjust the tightness of the latter.

- The height at which the stage rise can be limited using the locking lever situated at the base of the coarse-focusing knob on the left of the stand. To adjust the limit:
 - screw in the X 40 objective (or X 100 if one has been added).
 - lower the lever to unlock it by pushing it back;
 - gently raise the stage while looking at it sideways until the slide almost touches the lens (about 2 mm away);
 - raise the lever and lock it. The stage will not rise above the limit thus set. Repeat these operations to change the setting.
- Under the stage there is a light condenser that is moved using the knob on the right, an adjustable iris diaphragm and a filter holder designed to hold a (blue) filter for observing very transparent preparations.
- Adjust the lighting by seeking the best position of the condenser, then shut down the diaphragm to increase contrast and depth of field, but not too much so as to alter the resolution.
- Using the rotating stage: this is graduated to 2/10ths over 360° and is fitted with slide holders to hold histological or mineralogical slide preparations in place. A locking screw enables the stage to be locked in a chosen position.
- Centring the rotating stage: in order to align the rotational centre of the stage with the optical axis of the lens, the stage must be centred. If centring is not optimal it can be corrected using the two adjustment screws (on either side on the edge of the stage). Use a small screwdriver or suitable key.

- focus on the slide without polarizing the light, and mark out an easily identifiable point in the field of vision.
- position the slide so that this point lies in the centre of the eyepiece.
- rotate the stage; if the chosen point turns in a circle in the field of vision, the two screws must be adjusted to centre it correctly.
- gradually reduce the size of the circle in which the point moves until the latter remains fixed (see figure, below).
- the stage is now ready for observation of the slide with the objective in place.

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“Guide de pétrologie descriptive”, Ronan Hébert, Nathan Université, publishers.

4.1 Polarizer and analyser

The polarizer is placed above the light condenser; it is removable. It should be rotated 360° to lie at an angle of 90° to the analyser (see above: setting up, paragraph 6).

The analyser is incorporated in the head of the microscope. When the tab **(A)** is pushed in the analyser is not in place (light is not analysed); when the tab is pulled fully out, the analyser is in place (light is analysed).

5 Maintenance

5.1 Protection from dust

After use, cover the microscope with the dust-cover provided. If dust settles on the lenses, bow it off using a syringe or special duster. Lenses are delicate and can be scratched by dust.

5.2 Fingerprints

Do not touch the lenses. If they should become marked by fingers, remove the marks with a soft cloth or using a special lens cleanser and wipe dry with lens-cleaning tissue.

5.3 Damp


Moisture, especially in combination with changes in temperature, is particularly harmful to mounted lenses. Microscopes should be kept in a dry place in their polystyrene boxes with a desiccating sachet.

5.4 Dismantling

A microscope is a delicate precision instrument. If it has to be dismantled (for maintenance) it should be returned to the supplier or to sent to a qualified workshop.

The plate that closes it on the bottom should only be opened by an electrician, for the housing must remain electrically safe. Only dismantle the microscope when it has been unplugged from the mains.

5.5 Changing the bulb

1. Unplug the mains connection and remove the cable from the female socket at the back of the microscope.
2. Turn the microscope back to front and open the compartment housing the light by undoing the screw.
3. Ensure that the bulb has cooled.
4. Grasp the sides of the bulb and remove it.
5. Fit the new bulb taking care not to touch the lamp with your fingers.

6. Close the compartment and tighten the screw.
7. Return the microscope to its position.
8. Replace the flex on the microscope and plug into the mains socket.
9. Spare 6 V / 20 W light bulb.

5.6 Changing the fuse

Should the light fail to work even after changing the bulb, the fuse may need to be changed (spare fuse supplied).

This operation should be carried out by a qualified person (e.g. maintenance technician) and after the microscope has been disconnected from the mains electricity supply.

To access the fuse, unscrew the 4 screws in the plate on the bottom of the microscope.

Replace the defective fuse in the electronics chart with a new fuse.

Spare fuse: preferably use time-lag fuse, 250 mA, 250 V, 5 mm diameter, 20 mm long.


Only use approved fuses.

6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*