



Streamline Cell and General Separations

ALLEGRA® X-12 SERIES

BENCHTOP CENTRIFUGES

NEW



Streamline Separations • Save

Whether harvesting cells, separating blood components or other general separations — you can simplify your processes and save time with the Allegra® X-12 Series Benchtop Centrifuges. Offering a maximum capacity of 4 x 750 mL, these centrifuges specialize in streamlining cell and general separations.

ELIMINATE STEPS IN YOUR CULTURING PROCESS WITH EXCLUSIVE CELL CULTURE FLASK ADAPTERS*

By direct centrifugation of cell culture flasks, Allegra X-12 Series centrifuges eliminate the entire transfer step to and from conical tubes for cell separation:

- Saving Transfer Time and Labor
- Eliminating the use of 15 mL or 50 mL Conical Tubes
- Reducing a Potential Contamination Step

For more information on cell culture flask adapters, request AIB# A-1974A.



AUTOMATICALLY CORRECT IMBALANCE WITH ARIES™ SMART BALANCE SYSTEM


Exclusive ARIES 4 x 750 mL Self-Balancing Rotor:

- Automatically Corrects for up to 50 Grams of Imbalance in Opposing Loads – and Continues the Run Without Shutdown
- “Eye Balance” Samples Instead of Weighing – Saving Valuable Time



ENHANCE EASE-OF-USE WITH FEATURES NOT TYPICALLY FOUND IN BENCH CENTRIFUGES

- Numerical Key Pad for Easy & Quick Value Entry
- Control Temperature for up to 10 Minutes After the Rotor has Stopped, for Walk-Away Flexibility
- Automatic Initial “Soft Start” to Ensure Sample Integrity



FURTHER INCREASE CONVENIENCE WITH . . .

- 10 User-Defined Programs
- 10 Independent Acceleration/11 Independent Deceleration Profiles
- The Ability to Set by RPM, RCF and Toggle Between Display Modes
- Time Options: Set, Hold or Pulse (for Short Runs)
- Automatic Rotor Identification
- Audible End-of-Run Tone

OPTIMIZE THE USE OF POPULAR LABWARE FOR CELL AND GENERAL SEPARATIONS

4 x 750 mL Swinging Bucket Rotors (Standard & ARIES Self-Balancing) - 3,750 RPM; 3,270 x g

Adapters to Accommodate:

- Corning** 75cm² and 25cm² Canted-Neck Cell Culture Flasks
- 750 mL Flat-Bottom, 500 mL (Conical & Flat-Bottom) and 250 mL (Conical & Flat-Bottom) Bottles
- 15 mL (Conical & Round Tubes)
- 50 mL (Conical & Round Tubes)
- Blood Bags
- Vacutainer Tubes (3 mL, 5 mL, 7 mL and 10 mL)
- 1.5 mL Tubes
- And More




Save Time • Boost Convenience



BIOCONTAINMENT — VARIETY OF OPTIONS — EACH INDEPENDENTLY CERTIFIED BY CAMR, PORTON DOWN, UK

- Bio-Certified Covers for Round Buckets on 4 x 750 mL Standard and ARIES Rotors
 - Accommodate Specified Labware in Both Versions
- Exclusive Bio-Certified Aerosolve™ Canisters
 - Clear Plastic for Viewing from Top and Bottom
 - Compatible with 4 x 750 mL Standard and ARIES Rotors
- Bio-Certified 6 x 100 mL Fixed Angle Rotor



POPULAR 96-WELL LABWARE FOR CELLULAR ASSAYS AND QUICK MICROPLATE SPINS

Optional Microplate Carrier for Standard 4 x 750 mL Rotor (SX4750)

- 12 Standard Microplates, 4 Deep-Well/Square-Well Plates, 4 of 96-Tube Racks
- 3,750 RPM; 2,890 x g

CONTINUE THE SIMPLICITY IN DOWN STREAM APPLICATIONS

6 x 100 mL Fixed Angle Rotor — 10,200 RPM; 11,400 x g

- Further Isolate Biomolecules
- Clearing of Media, Cell Debris, Lysate, Tissue Homogenates
 - 100 mL, 94 mL, 85 mL Bottles
 - 50 mL & 15 mL (Conical & Round) Tubes
 - Centrifugal Filter Devices
 - 1.5 mL Tubes
 - And More



AVAILABLE IN BOTH REFRIGERATED (ALLEGRA X-12R) AND CONSTANT CONTROLLED TEMPERATURE (ALLEGRA X-12) VERSIONS.

The Allegra X-12R Refrigerated Centrifuge allows selection of temperature within the range of -10° C to 40° C.

The Allegra X-12 Constant Controlled Temperature Centrifuge is preset to 20° C and automatically controls temperature within the range of 18° C to 22° C. Versus a traditional "non-refrigerated" centrifuge, you can be assured your sample will run no lower than 18° C and no higher than 22° C.

SPECIFICATIONS

	ALLEGRA® X-12R	ALLEGRA® X-12
Refrigeration	Yes, CFC-Free	Constant Controlled Temperature
Maximum Speed	10,200 rpm (Fixed-angle)	10,200 rpm (Fixed-angle)
Maximum g-Force	11,400 x g (Fixed-angle)	11,400 x g (Fixed-angle)
Maximum Capacity	4 x 750 mL	4 x 750 mL
Run Time	Timed (Up to 99 hours, 59 min) Hold, Pulse (Short Run)	Timed (Up to 99 hours, 59 min) Hold, Pulse (Short Run)
Drive	Brushless Induction	Brushless Induction
Noise at Max Speed	≤68 dBA	≤68 dBA
User-Defined Programs	10	10
Accel/Decel Profiles	10/11	10/11
Temperature Range	-10°C to +40°C (Selectable)	+18°C to +22°C (Preset)
Dimensions(cm)	76.2W x 62.2D x 34.3H	76.2W x 62.2D x 34.3H
Dimensions(in)	30W x 24.5D x 13.5H	30W x 24.5D x 13.5H
Weight without Rotor	121 kg (267 lbs)	121 kg (267 lbs)

ORDERING INFORMATION

PART No.	DESCRIPTION	MAX SPEED(rpm)	MAX RCF, g
		SWINGING BUCKET/FIXED ANGLE	SWINGING BUCKET/FIXED ANGLE
392302	Allegra X-12R Centrifuge, 208V, 60-Hz	3,750 / 10,200	3,270 / 11,400
392304	Allegra X-12R Centrifuge, 230V, 50-Hz	3,750 / 10,200	3,270 / 11,400
392303	Allegra X-12R Centrifuge, 200V, 50/60-Hz	3,750 / 10,200	3,270 / 11,400
392472	Allegra X-12 Centrifuge, 208V, 60-Hz	3,750 / 10,200	3,270 / 11,400
392474	Allegra X-12 Centrifuge, 230V, 50-Hz	3,750 / 10,200	3,270 / 11,400
392473	Allegra X-12 Centrifuge, 200V, 50/60-Hz	3,750 / 10,200	3,270 / 11,400
SWINGING BUCKET ROTORS			
369702	4x750 mL Swinging Bucket Rotor Assembly w/o Bucket Covers (SX4750)	3,750	3,270
369704	4x750 mL ARIES Self-Balancing Swinging Bucket Rotor Assembly w/o Bucket Covers (SX4750A)	3,750	3,270
392805	Bio-Certified† Covers for SX4750 and SX4750A Round Buckets (Set of 2)		
392806	Microplate Carriers for SX4750 Rotor (SX4750μ)(Set of 2) Includes Carriages (PN392873)	3,750	2,890
392807	Covers for SX4750μ Microplate Carriers (Set of 2)		
392873	Replacement Carriages for SX4750μ Microplate Carriers, Includes Rubber Pads (Set of 2)		
FIXED ANGLE ROTOR			
369735	Bio-Certified† 6X100mL Fixed Angle Rotor (FX6100)	10,200	11,400

PART No.	VOLUME	DIAMETER(mm)	COLOR	BUCKET CAPACITY	ROTOR CAPACITY	PACK OF
ADAPTERS FOR DIRECT CENTRIFUGATION OF CELL CULTURE FLASKS (SX4750 & SX4750A)						
369292	Corning 75-cm ² Canted-Neck Cell Culture Flasks	NA	Orange	1	4	2
369295	Corning 25-cm ² Canted-Neck Cell Culture Flasks	NA	Green	2	8	2
MODULAR DISK ADAPTERS FOR TUBES (SX4750 & SX4750A)						
359469 & 359148	3 mL & 5 mL	10 & 12	Blue	37	148	2 & 4
343108	Tube Decanter for PN 359469 & 359148	NA	NA	NA	NA	Each
354511	1.5 mL Adapter Plate for PN 359469 & 359148	NA	NA	NA	NA	Each
359478 & 359157	3 mL & 5 mL	13	Tan	30	120	2 & 4
359470 & 359149	7 & 10 mL	14	Orange	24	96	2 & 4
359471 & 359150	12 mL	16	Purple	19	76	2 & 4
359472 & 359151	15 mL Conical	18	Green	14	56	2 & 4
359473 & 359152	15 mL & 20 mL Round	18	Green	14	56	2 & 4
359475 & 359154	30 mL Round & 50 mL Conical	30	Light Green	4	16	2 & 4
359474 & 359153	50 mL Round	29	Yellow	7	28	2 & 4
359476 & 359155	50 mL Round	35	Dark Blue	4	16	2 & 4
359477 & 359156	100 mL	44	Brown	2	8	2 & 4
ADAPTERS FOR BOTTLES & BLOOD BAGS (SX4750 & SX4750A)						
356983	230 mL Conical (Insert for Use with 349946)	62	Orange	1	4	Each
349946	250 mL Flat Bottom	62	Yellow	1	4	Each
349849	250 mL Conical	62	Warm Red	1	4	Each
349945	500 mL Flat Bottom	70	Tan	1	4	Each
349945	500 mL Conical (with Corning Adapter)	70	Tan	1	4	Each
349846	750 mL Flat Bottom	96	Blue	1	4	Each
356856	Single or Double Pack Blood Bag	90	Yellow	2	8	Each
356857	Triple or Quad Pack Blood Bag	97	Orange	1	4	Each
BIO-CERTIFIED† AEROSOLVE CANISTERS AND TUBE RACKS (SX4750 & SX4750A)						
359481 & 359232	Bio-Certified Aerosolve Canisters	NA	NA	NA	NA	2 & 4
354495	1.5 mL	11	White	24	96	Each
359482 & 359160	3 & 5 mL	12	Blue	24	96	2 & 4
359489 & 358993	5 mL	12	Tan	24	96	2 & 4
359483 & 359161	10 mL	14	Orange	18	72	2 & 4
359484 & 359162	12 mL, 3 & 5 mL	16,12	Purple	12,6	48,24	2 & 4
344517	15 mL Round	14	White	10	40	Each
359485 & 359163	15 & 20 mL Round, 3 & 5 mL	18,12	Green	12,6	48,24	2 & 4
359487 & 358991	15 mL Conical, 3 & 5 mL Round	17,12	Light Green	6,6	24,24	2 & 4
359486 & 359164	50 mL Round, 3 & 5 mL	29,12	Yellow	4,4	16,16	2 & 4
359488 & 358992	50 mL Conical, 3 & 5 mL Round	30,12	Lime Green	4,4	16,16	2 & 4
356985	230 mL Conical	62	Orange	1	4	2
ADAPTERS FOR TUBES (FX6100)						
344497	1.5/2.0 mL (must be inserted into PN 392830)	NA	White	NA	36	6
392824	10 mL Round	16	White	NA	6	6
392823	15 mL Round	18	White	NA	6	6
392270	15 mL Conical	17	Beige	NA	6	6
392822	30 mL Round	25.7/25	White	NA	6	6
392821	35 mL Round	29	White	NA	6	6
392830	50 mL Round	29	Blue	NA	6	6
392268	50 mL Conical	28.5	Beige	NA	6	6

Accommodates 100 mL, 94 mL, 85 mL, 80 mL and 70 mL tubes without adapters, directly in rotor cavity.

* Patent Pending

** Corning is a trademark of Corning Incorporated.

†Certified by CAMR, Porton Down, UK for Biocontainment.



Developing innovative solutions in genetic analysis, drug discovery, and instrument systems.

Innovate Automate
SIMPLIFY

Beckman Coulter, Inc. • 4300 N. Harbor Boulevard, Box 3100 • Fullerton, California 92834-3100

Sales & Service: 1-800-742-2345 • Telex: 678413 • Fax: 1-800-643-4366 • www.beckmancoulter.com

Worldwide Biomedical Research Division Offices:

Australia (61) 2 9844-6000 Canada (905) 819-1234 China (86) 10 6515 6028 Eastern Europe, Middle East, North Africa (41) 22 994 07 07

France 01 49 90 90 00 Germany 49 21 513335 Hong Kong (852) 2814 7431/2814 0481 Italy 02-953921 Japan 03-5404-8359

Mexico (55)-560-57770 Netherlands 0297-230630 Singapore (65) 6339 3633 South Africa, Sub-Saharan Africa (27) 11-805-2014/5 Spain 91 3836080

Sweden 08-564 85 900 Switzerland 0800 850 810 Taiwan (886) 2 2378 3456 Turkey 90 216 309 1900 U.K. 01494 441181 U.S.A. 1-800-742-2345