

BENEFITS

- Fast, efficient cryogenic grinding at -196 °C
- No bothersome manual cooling thanks to integrated cooling system with Autofill
- High final fineness due to grinding by impact and friction
- Highly reproducible grinding results
- Programmable cooling and grinding cycles
- Digital parameter setting
- Simple operation, easy cleaning
- Feed size: up to 8 mm
- Final fineness: down to approx. 5 μm
- With 1 grinding station for jars up to 50 ml www.retsch.com/cryomill

IDEAL FOR:

- PLASTICS
- TEMPERATURE-SENSITIVE SAMPLES WITH VOLATILE
 - COMPONENTS

Referen

Solutions in Milling & Sieving

PARTICULARLY SAFE:

- AUTOFILL SYSTEM FOR LIQUID NITROGEN
 - NO-LOSS SCREW-TOP GRINDING JARS UP TO 50 ML

CRYOMILL



1	Performance data	CryoMill		
	Field of application	size reduction, mixing, homogenization, cell disruption		
	Feed material	hard, medium-hard, soft, brittle, elastic, fibrous		
	Feed size*	up to 8 mm		
	Final fineness*	approx. 5 μm		
	Batch/Sample volume	max. 1 x 20 ml / 4 x 2 ml		
	Typical mean grinding time	pre-cooling: 10 minutes, grinding: 1 minute		
	Possible applications			
	Cryogenic grinding	yes		
	Wet grinding	yes		
	Dry grinding	yes		
	Cell disruption	yes		
	Suitable grinding jars	5, 25, 35, 50 ml		
	No. of grinding stations	1		
	Digital preselection of vibrational freque	ency 3 - 25 Hz (180 - 1500 min ⁻¹)		
	Digital preselection of grinding time	10 s - 99 min		
	Memory for parameter combinations	1		
4				
	Technical data			
	Power consumption	160 W		
	$W \times H \times D$	385 x 370 x 570 mm		
	Net weight	approx. 46 kg		

*depending on feed material and instrument configuration/settings

Order data

ı	CryoMill				Item No.		
	CryoMill						
	(Please order Autofill with LN ₂ container, grinding jars and balls separately)						
7	CryoMill for 100–240 V	V, 50/60 Hz			20.748.0001		
	Grinding jars	5 ml*	25 ml	35 ml	50 ml		
	Stainless steel	01.462.0290	01.462.0289	01.462.0288	01.462.0284		
	*to be used with adapter 02.706.0272						
	Accessories						
ı	Autofill with ${\rm LN_2}$ container, 10 l	02.480.0001					
ŀ	Adapter for use of 2/4 grinding	02.706.0272					

For additional accessories as grinding balls etc. please refer to RETSCH's brochure "Ball Mills".

With integrated cooling system

The CryoMill has been specially designed for cryogenic grinding. The grinding jar is continually cooled with liquid nitrogen from the integrated cooling system before and during the grinding process. Thus the sample is embrittled and volatile components are preserved. The liquid nitrogen circulates through the system and is continually replenished from an autofill system in the exact amount which is required to keep the temperature at -196 °C. The automatic cooling system guarantees that the grinding process is not started before the sample is thoroughly cooled. This results in reduced consumption and guarantees reproducible grinding results.

The grinding parameters can be stored which helps to simplify routine tasks. LEDs in the display indicate the current state of operation, e.g. cooling or grinding.

The size reduction principle is the same as that of the MM 400. With a **vibrational frequency of 25 Hz** the CryoMill grinds most materials very effectively in a few minutes. The sample is mainly ground by impact but also by friction which allows for **substantially finer grind sizes** compared to other cryogenic mills.

The CryoMill is equipped with one grinding station for grinding jar volumes of 25 ml, 35 ml and 50 ml. Another option is the use of an adapter which holds up to four 5 ml grinding jars.

The mill can also be operated without cooling which makes it suitable for a vast range of applications.



Retsch GmbH

Rheinische Straße 36 42781 Haan, Germany

 Telephone
 +49 2129/5561-0

 Telefax
 +49 2129/8702

 E-mail
 info@retsch.com

 Internet
 www.retsch.com

a VERDER company

RETSCH – Your specialist for sample preparation offers you a comprehensive range of equipment. Please request information on our crushers, mills, sieve shakers, sample dividers, feeders as well as cleaning and drying machines.