


PuraLev®

Life Science Integrated Pump Series



PuraLev® i30SU (Single-Use)

| | |
|----------------|-------------------|
| 1.0 bar | (14.5 psi) |
| 7.7 liters/min | (2.0 gallons/min) |

Low Shear Design - High Cell Viability




Figure 1: Schematic of the main elements of the MagLev centrifugal pump




Figure 2: Integrated MagLev pump driver with pump head




Figure 3: Single-use pump head concept

INTRODUCTION

Levitronix® has developed a revolutionary pump that has no bearings to wear out or seals to break. Based on the principles of magnetic levitation, the pump's impeller is suspended, contact-free, inside a sealed casing and is driven by the magnetic field of the motor (Figure 1). The impeller and casing are both fabricated from biocompatible (FDA, USP-VI, BSE/TSE and Animal free) gamma sterilizable polypropylene (PP) and together they make up the disposable pump head. A simple and intuitive exchange of the single use pump head is achieved with a bayonet socket type mounting procedure (see Figure 6). Flow rate or pressure is precisely controlled by electronically regulating the rotor speed, which eliminates any pulsation. With the lack of mechanical bearings plus the self-contained pump head design, the risk of contamination is drastically reduced. The absence of narrow gaps between the impeller and pump casing, plus the low-shear pump design allows the gentle pumping of sensitive liquids. The pump head can be easily inserted and removed with an intuitive bayonet socket. The controller and the motor are integrated into the driver housing (see Figure 2), hence cabling effort is reduced.

SYSTEM BENEFITS

- Low shear-forces.
- Reduced risk of contamination due to the self-contained design with magnetic bearings.
- No particle generation.
- No over-pressure situations (compared to roller pumps).
- No narrow gaps between the impeller and pump casing where bacteria could be entrapped.
- Pump head is gamma sterilizable.
- Biocompatibility of wet materials: FDA, USP-VI, Animal/BSE/TSE free.
- Bayonet socket design for easy and intuitive exchange of disposable pump head (see Figure 6).
- Small size.
- Dry running capability.
- Proven technology in the medical (disposable blood pumps) and semiconductor (high-purity pumps) industries.
- High flow capability with compact design.
- Pulsation free.

APPLICATIONS

- Pumping of shear-sensitive liquids and cells.
- Bioprocessing (for example perfusion).
- Recirculation and transfer applications in bioreactors.
- Filtration.

SYSTEM CONFIGURATION – “STAND-ALONE”

Figure 7 and Figure 11 illustrate a “Plug and Play” stand-alone system with integrated user panel and buttons to set the speed manually.

The driver also contains a PLC interface for remote speed control by analog and digital signals.

Various accessories are available like a desktop power supply with relevant power cable and signal cables to connect to the PLC.

SYSTEM CONFIGURATION – “EASYCONNECT”

The “EasyConnect” models (see Figure 8 and Figure 13) with according cable accessories are designed to realize various interface configurations with minimal setup effort.

Two Fieldbus connectors (IN and OUT) allow to setup arrays of multiple pumps. Therefore, serial pumping configurations as shown in Figure 9 can be realized.

The PLC interface allows not only remote control by analog/digital signals but also connections of external sensors hence enabling for example a precise flow or pressure control (see notes below).

The Fieldbus interface allows remote control over a PC, a User Panel or other devices with Modbus protocol.

SYSTEM CONFIGURATION – “OEM”

The “OEM” models are designed for a compact integration with one integrated driver cable containing all available interface signals (see Figure 10 and Figure 15). Basically all configurations of the “EasyConnect” models are possible allowing the users with integration capabilities to adapt the cable to their needs.

PROCESS CONTROL WITH FEEDBACK SENSORS

Together with an external sensor, process parameters like flow or pressure can be controlled or monitored as shown in Figure 10.

Precise ultrapure flow control systems can be realized with the pump system in combination with LEVIFLOW® flowmeters. Levitronix® provides either turnkey solutions for closed loop flow control or helps to design your own flow control system. Experience has been gained a wide range of applications.

The versatility of Levitronix® flow control systems goes far beyond the capabilities of simple flow controllers. In addition to the flow control function, the Levitronix® control firmware comes with several condition monitoring features to monitor the integrity of the fluid circuit. Levitronix® flow control systems can generate alarms for preventive filter exchange, no-flow conditions or line clogging. Dynamic Condition Trending (DCT) enables failure prediction and scheduling of preventive maintenance.




Figure 4: Pressure/flow curves for aqueous liquids (similar to water)




Figure 5: Pump system models




Figure 6: Intuitive 3-step pump head mounting procedure with bayonet type socket (PHS-i30.1)




Figure 7: Standard "Stand-Alone" system configuration with main accessories
(See section "Order Information" for details to numbered components and other options)




Figure 8: Standard "EasyConnect" system configuration with main accessories
(See section "Order Information" for details to numbered components and other options)




Figure 9: Serial pumping configuration with "EasyConnect" models
(See section "Order Information" for details to numbered components and other options)




Figure 10: Standard "OEM" system configuration
(See section "Order Information" for details to numbered components and other options)



| Interface | PIN Name | Description | Standard Designation | Hardware Specification |
|---------------------|----------|------------------------------|--------------------------|---|
| Power Supply | P+ | + 24 VDC | Supply | Voltage: 24 VDC |
| | P- | Power Input Ground / Earth | | Power: 35 W |
| | NC | Not connected. | -- | -- |
| PLC 6 | Ain | Analog Input (Current Input) | Remote Speed | Analog current input; 4 – 20 mA (450 Ohm shunt input, no galvanic isolation) |
| | Ain_GND | Analog In. GND | -- | Reference for Ain |
| | Dout | Digital Output 1 | Status | Open drain, max. 24V, 100mA Reference ground is GND |
| | GND | Analog Ground | -- | Reference for Dout |
| | Din1 | Digital Input 1 | Enable (Reset) | Galvanic separation with optocoupler 2.2 kΩ input resistance, 5-24V for active input |
| | Din_COM | Com. Digi. Input | -- | Reference for digital input. |
| Display and Buttons | -- | Display | Speed and Status Display | -- |
| | -- | Up/Down | Setting speed | -- |
| | -- | On/Off | Enable/Disable | -- |

Figure 11: Interface specifications of standard "Stand-Alone" model






Figure 12: Basic dimensions and description of standard "Stand-Alone" model

MODEL DESCRIPTION – EASYCONNECT




| Connector | PIN Name | Description | Standard Designation | Hardware Specification | |
|--------------|-------------|--------------------------------|----------------------|--|--|
| Power Supply | P+ | + 24 VDC | Supply | Voltage: 24 VDC | |
| | P- | Ground / Earth | | Power: 35 W | |
| | NC | Not connected. | -- | -- | |
| PLC 12 | Dout1 | Digital Output 1 | Status | Open drain, max. 24V, 100mA | |
| | Dout2 | Digital Output 2 | Error | Reference ground is GND | |
| | Din1 | Digital Input 1 | Enable (Reset) | Galvanic separation with optocoupler | |
| | Din2 | Digital Input 2 | Process Mode | 2.2 kΩ input resistance, 5-24V for active input | |
| | Din_COM | Com. Digi. Input | -- | Reference for digital input. | |
| | Ain1 | Analog Input 1 (Current Input) | Actual Process Value | Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation) | |
| | Ain2 | Analog Input 2 (Voltage Input) | Reference Value | Analog voltage input: 0 – 10V (7.9 kΩ, no galvanic isolation) | |
| | Ain_GND | Analog In. GND | -- | Reference for Ain1 and Ain2 | |
| | Aout1 | Analog Output (Voltage Output) | Actual Speed | 0 – 10V (no galvanic isolation) GND is reference | |
| | GND | Analog Ground | -- | Reference for Aout1, Dout1, Dout2 and Pout | |
| | Pout | Output +24VDC | Supply Output | For supply of external devices (e.g. sensors). (Current 200mA together with Pout o Fieldbus OUT) | |
| | NC | Not connected. | -- | -- | |
| Fieldbus OUT | GND | Ground | -- | Reference for Pout. | |
| | Pout | Output +24VDC | Supply Output | For supply of external devices (user panels) (Current 200mA together with Pout of PLC 12) | |
| | RS485+ | RS485 + | Field Bus | Modbus protocol | |
| | RS485- | RS485 - | | | |
| | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. | |
| | Internal | Internal Bus | Do not connect | | |
| | GND | Ground | -- | Connected to PLC12 GND and reference for supply | |
| | NC | Not connected. | -- | -- | |
| | Fieldbus IN | RS485+ | RS485 + | Field Bus | Modbus protocol |
| | | RS485- | RS485 - | | |
| | | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. |
| | | Internal | Internal Bus | Do not connect | |


Figure 13: Interface specifications of standard "EasyConnect" model



Pump Head Inlet and Outlet Fitting:
Barb 3/8" for tubing with typical ID = 1/4"




Single-Use Pump Head DCP-30.2



Note: Pump Head Socket can be mounted and rotated in steps of 45° to change pump outlet orientation.

Figure 14: Basic dimensions and description of standard "EasyConnect" model



| Wire Name | Description | Standard Designation | Hardware Specification |
|-----------|--------------------------------|----------------------|---|
| P+ | + 24 VDC | Supply | Voltage: 24 VDC |
| P- | Power Input Ground / Earth | | P- to be connected to earth |
| Ain1 | Analog Input 1 (Current Input) | Actual Process Value | Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation) |
| Ain2 | Analog Input 2 (Voltage Input) | Reference Value | Analog voltage input: 0 – 10V (7.9 kOhm, no galvanic isolation) |
| Ain_GND | Analog Input Ground | -- | Reference for Ain1 and Ain2 |
| Din1 | Digital Input 1 | Enable (Reset) | Galvanic separation with optocoupler 2.2 kΩ input resistance, 5-24V for active input |
| Din2 | Digital Input 2 | Process Mode | |
| Din_COM | Common Digital Input | -- | -- |
| Aout1 | Analog Output (Voltage Output) | Actual Speed | 0 – 10V (no galvanic isolation) GND is reference |
| Dout1 | Digital Output 1 | Status | Open drain, max. 24V, 100mA |
| Dout2 | Digital Output 2 | Error | Reference ground is GND |
| GND | Analog Ground | -- | Reference for Aout1, Dout1 and Dout2 |
| RS485+ | RS485 + | Field Bus | Modbus protocol |
| RS485- | RS485 - | | |
| Internal | Internal Bus | Do not connect | For internal usage. |
| Internal | Internal Bus | Do not connect | For internal usage. |
| Shield | Shielding | Shielding | To be connected to earth (see wire No. 2, P-) |

Figure 15: Interface specifications of standard "OEM" model

Note 1: Power supply wire cross-section is 1.5 mm² and for signal wires 0.14 mm²

Note 2: For more detailed description of interfaces consult user manual






Figure 16: Basic dimensions and description of standard "OEM" model


7 IPS Cable Power ICP-1.1 (for EasyConnect and Stand-Alone drivers)



8a IPS Cable Signal ICS-1.1 (for EasyConnect drivers)



8b IPS Cable Signal ICS-1.2 (for EasyConnect and Stand-Alone drivers)



9 IPS Cable Signal ICS-2.1 (for EasyConnect drivers)




Figure 17: Basic dimensions and specifications of standard cables

ORDER INFORMATION

| System Name | Article # | Pump Head Socket | Driver | Note |
|-------------|-----------|------------------|----------------------------------|--|
| PLD-i30SU.1 | 100-90955 | PHS-i30.1 | IPD-30.1-50-02 | OEM - Driver, 5 m PVC cable with open wires, pump head socket |
| PLD-i30SU.2 | 100-91025 | PHS-i30.1 | IPD-30.3-02 (MBP-i30.1 included) | EasyConnect - Driver with interface connectors, pump head socket. |
| PLD-i30SU.3 | 100-91026 | PHS-i30.1 | IPD-30.5-02 (MBP-i30.1 included) | Stand-Alone - Driver with integrated user panel, pump head socket. |

Table 1: Standard driver system configurations

| Pos. | Component | Article Name | Article # | Characteristics | Value / Feature |
|----------|---|--|------------------------|---|---|
| 1a 1b | Single-Use (SU) Pump Head | DCP-30.2 (Barb) DCP-30.1 (Triclamp) | 100-90968 100-90959 | Material Impeller and Pump Housing Housing Sealing In-/Outlet Fittings Max. Flow Max. Diff.-Pressure Max. Viscosity Max. Liquid Temp. Wet Pump Volume/Surface Sterilization Methods | Polypropylene (FDA, USP Class VI, BSE/TSE/Animal free) Infrared welding Barb 3/8" or Triclamp 3/8" for tubing with typical ID = 1/4" 7.7 liters/min / 2.0 gallons/min 1.0 bar / 14.5 psi 40 cP 60 °C / 140 °F 7.7 ml / 55.9 cm ³ Gamma radiation up to 40kGy |
| 1c 1d | Irradiated SU Pump Head | DCP-30.2-G25 (Barb) DCP-30.1-G25 (Triclamp) | 100-91071 100-91170 | Applied Gamma Dosage | ≥ 25 kGy |
| 1e | Irradiated SU Pump Head with Sterile Fittings | DCP-30.2-SF1-G25 | 100-91234 | Pump Type (A) / Tubing (C) Sterile Fittings (B) Fitting Compatibility Applied Gamma Dosage | DCP-30.2 (Barb) / Silicone AseptiQuik® S from CPC® with part # AQS17006 ³ Various sizes and types including autoclavable versions available at CPC® ≥ 25 kGy |
| 2a | Integrated Pump Driver ("OEM Model") | IPD-30.1-50-02 | 100-10088 | Voltage, Power Housing Cable Interfaces Standard Firmware | 24 VDC ±10%, 35 W Epoxy coated Aluminum, PP for bottom lid, IP65 ¹ PVC jacket, open wires, cable length 5 m PLC and RS485 with Modbus protocol (see Figure 15 for details) H2.48 |
| 2b | Integrated Pump Driver ("EasyConnect" Model) | IPD-30.3-02 (MBP-i30.1 included) | 100-10097 | Housing Interfaces Standard Firmware | Epoxy coated Aluminum, PP for bottom lid, IP65 2x Fieldbus RS485 with Modbus protocol, PLC and power supply H2.48 ² |
| 2c | Integrated Pump Driver ("Stand-Alone" Model) | IPD-30.5-02 (MBP-i30.1 included) | 100-10098 | Housing Interfaces Standard Firmware | Epoxy coated Aluminum, PP for bottom lid, IP65 User panel with 3 user buttons, PLC and power supply H2.25 |
| 3 | Pump Head Socket | PHS-i30.1 | 100-90947 | Mounting Type Material Assembly Screws | Bayonet type with locking pin Anodized Aluminum 4 pcs M3 x 6 mm (Stainless Steel, INOX A4) |

Table 2: Specification of standard components

Note 1: Designed and tested for IP67. Note 2: Special firmware for serial pumping (see Figure 9) as one unit available on request. Note 3: CPC® and AseptiQuik® are registered marks of the Colder Product Company.

| Pos. | Component | Article Name | Article # | Characteristics | Value / Feature |
|------|---|--|---|---|---|
| 4 | Mounting Base Plate | MBP-i30.1 | 190-10313 | Material / Mounting Screws | PP + 30% GF / 2 pieces, stainless steel FEP coated, M3 x 10 mm |
| 5a | AC/DC Power Supply | TPC 055-124 (Traco) | 100-40014 | Voltage Output / Input Basic Dimensions Certification or Standards | 24 VDC with 55 W / 85 – 264 VAC, 47-63 Hz 45 x 90 x 96.5 mm (mountable on DIN rail 35 mm) UL, CSA, CB, Semi F47 |
| 5b | Desktop AC/DC Power Supply | AC/DC Power Supply VEC50US24 HR30 | 100-40015 | Voltage Output / Input Basic Dimensions Safety Approvals Note | 24VDC, 50W / 90 – 264 VAC, 47-63 Hz 116 x 52 x 31 mm IEC60950-1, EN60950-1, UL/cUL60950-1 Connector for direct connection to power supply of driver with cable length 1.2m. |
| 5c | AC Mains Cables (for Desktop power supply 5b) | AMC-1.1 (2m) AMC-1.2 (2.5m) AMC-1.3 (2.5m) AMC-1.4 (2.5m) AMC-1.5 (2.5m) | 190-10331 190-10332 190-10333 190-10334 190-10335 | Approvals and Country Approvals and Country Approvals and Country Approvals and Country Approvals and Country Cable Specifications | UL, cUL, US, Canada CB, Germany, Denmark, Norway, Finland, Belgium, Netherland, Sweden, Austria PSE, Japan Switzerland CE, United Kingdom Black color, ROHS |
| 6 | USB to RS485 Adaptor-TR Isolated | YN-485I-TR | 100-30392 | Structure/Design Purpose | USB connector (A) with termination resistor and cable (2m) with connector pair (B and C) for external RS485 wire connection. Magnetically isolated. Cable length is 2m. Included is a USB space saver cable (D). Communication over fieldbus of driver with PC |
| 7 | IPS Cable Power 3 Wires | ICP-1.1-50 (5 m) | 190-10342 | Cable Material / Wires Connection In / Connection Out Main Purpose | PVC jacket / 3x 0.5 mm ² (only 2 wires used, 1 is cut) Open wires / Circular Hirose type to driver Connection of power supply to "Stand-Alone" and "EasyConnect" drivers |
| 8a | IPS Cable Signal 6 Wires | ICS-1.1-01 (0.1 m) ICS-1.1-10 (1 m) ICS-1.1-30 (3 m) | 190-10343 190-10344 190-10345 | Cable Material / Wires Connection In / Connection Out Main Purpose | PVC jacket / 6x 0.08 mm ² and shielding Circular Hirose type / Circular Hirose type Fieldbus connection between "EasyConnect" drivers (e.g. multi-pump arrays) |
| 8b | IPS Cable Signal 6 Wires | ICS-1.2-50 (5 m) | 190-10346 | Cable Material / Wires Connection In / Connection Out Main Purpose | PVC jacket / 6x 0.08 mm ² and shielding Connector with screw type plug for open wire connection / Circular Hirose type Fieldbus connection to "EasyConnect" drivers and to PLC of "Stand-Alone" drivers. |
| 9 | IPS Cable Signal 12 Wires | ICS-2.1-50 (5 m) | 190-10347 | Cable Material / Wires Connection In / Connection Out Main Purpose | PVC jacket / 12x 0.14 mm ² and shielding Connector with screw type plug for open wire connection / Circular Hirose type General connection to PLC of "EasyConnect" drivers. |
| 10 | Fieldbus Termination Connector | FTC-1.1 | 190-10348 | Materials Main Purpose | PPS for connector housing and FPM for sealing. Termination of fieldbus. |
| 11 | User Panel | LUI-B.1-01 | 100-30448 | Interface / Housing Rating Standard Firmware | RS485 / IP65 A3.00 |

Table 3: Specification of accessories

Note 1: Designed and tested for IP67.



Figure 18: Pump systems with standard components



Figure 19: General standard accessories



Figure 20: Standard cables

Levitronix® is the world-wide leader in magnetically levitated bearingless motor technology. Levitronix® was the first company to introduce bearingless motor technology to the Semiconductor, Medical and Life Science markets. The company is ISO 9001 certified. Production and quality control facilities are located in Switzerland. In addition, Levitronix® is committed to bring other highly innovative products like the LEVIFLOW® flowmeter series to the market.



Headquarter and European Contact

Levitronix GmbH
Technoparkstr. 1
CH-8005 Zurich
Switzerland

Phone: +41 44 445 19 13
Fax: +41 44 445 19 14
E-Mail: salesEurope@levitronix.com

US Contact

Levitronix Technologies Inc.
10 Speen Street, Suite 102
Framingham, Massachusetts 01701
USA

Phone: +1 508 861 3800
Fax: +1 508 861 3837
E-Mail: salesUS@levitronix.com

Japan Contact

Levitronix Japan K.K.
Wing Eight 5floor, 4-16-4
Asakusabashi, Taito-ku
Tokyo, 111-0053 Japan

Phone: +81 3 5823 4193
Fax: +81 3 5823 4195
E-Mail: salesJapan@levitronix.com

Taiwan Contact

Levitronix Taiwan
5F, No. 251, Dong Sec. 1,
Guangming 6th Rd., Chu Pei City,
Hsin-Chu 302, Taiwan, R.O.C.

Phone: +886 3 657 6209
Fax: +886 988 321472
E-Mail: salesAsia@levitronix.com