Microscopie

Microscopy

**Ref :
571 168**

Français – p 1

English – p 7

Version : 6009

Microscope ME 1012

ME 1012 microscope

Microscopie
Microscope ME 1012
Ref :
571 168

JEULIN

FRANÇAIS

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JEULIN

1 General characteristics

The ME 1012 microscope is designed for dry observation of preparations in transmitted light on a light background.

- Stand with a modern, sturdy, stable and ergonomic design
- Low position of focusing controls allowing the user to keep his hands resting on the table
- Eyepiece holder head tilted at 45°: comfort during observation, stability of preparations in liquid medium on the horizontal stage
- Head revolving around a vertical axis
- Revolving turret for 3 objective lenses
- Optical length 160 mm
- Rectangular mounting plate 110 x 120 mm equipped with 2 spring clamps and lever, fixed by a Phillips head screw
- Focusing by tooth rack and micrometer screw for vertical displacement of the stage
- Fixed condenser, F number 0.65 fixed to the stage, equipped with an iris aperture, a filter holder and a removable blue filter
- Lighting device built into the base, equipped with a 12 volts 10 watts halogen lamp, fitted with a permanent blue filter, a condensing lens and a switch
- Power supplied by a 12 volts transfer built into the mains socket
- Magnification of 40 x to 400 x
- Weight of the microscope: 3.3 kg.

2 Optical characteristics

Wide field eyepiece and index marker	10 x		
Achromatic objectives			
Magnification	4 x	10 x	40 x
F number	0.1	0.25	0.65
Mount	fixed	fixed	telescopic
Working distance	16.00 mm	7.31 mm	0.45 mm
Field of view	4.00 mm	1.65 mm	0.40 mm

3 Commissioning

Take the microscope out of the polystyrene cabinet.

Install the three objectives on the turret.

Reserve the removable blue filter for special use: observing particularly transparent objects at 40 magnifications.

Cover the device with its dust cover if it is not used immediately.

Otherwise, connect the 12 volts plug at the back of the microscope and connect the 220 volts transformer plug to the mains: the microscope is ready for operation.

4 Use and adjustment

4.1 Clasp bar and base

The clasp bar with a large section provides a tilted grasping area that is very convenient for moving the microscope: it is strongly advised not to use other parts, which include precision mechanisms, for this purpose.

The stable base supports the column and the lighting device with its switch. The lower plate that can be closed must be opened only by an electrician, as the box must remain electrically secure.

4.2 Lighting

Toggle the switch: the lamp comes on. It includes a 12V / 10W halogen lamp with a dichroic mirror, a blue correction filter (eliminating the mainly red and yellow colours of the incandescent lamp) and a converging lens on the axis of the condenser.

To change the lamp:

1. Disconnect the power supply, either using the 12 volts plug located at the back of the microscope, or by the 220 volts mains plug
2. Completely unscrew the black lamp body
3. Make sure that the lamp has cooled down
4. Avoid touching the bulb
5. Hold the outer edges of the mirror and pull to remove the lamp
6. Replace and screw back on the lamp body before connecting the power supply.

Spare lamp 12V / 10 W part no.: 574202.

4.3 Observation head

The head is tilted at 45° to obtain better comfort for observation. It can pivot by 360° around a vertical axis, which allows the user to make observations in all directions without having to move the base.

4.4 Eyepiece and objectives

The eyepiece with a wide field, magnification 10, has a permanent directed index marker at its centre. By rotating the body of the eyepiece, this index marker can be placed in the position that is most favourable for observation. The eyepiece is fixed with a screw, which prevents it from accidentally falling off during transport.

The achromatic objectives, installed correctly on the turret, are of the parfocal type: it is therefore not necessary to significantly adjust the focus when changing an objective.

The 40-x objective is of the telescopic type: its spring mount prevents damage and breaking of the preparation in case of accidental contact.

An adjustable stop helps, moreover, prevent this contact (see "focusing").

4.5 Stage

The rectangular 110 x 120 stage is equipped with two clamps installed on a spring fixed by a screw preventing them from falling off and getting lost. Each of them tips like a lever simply by pressing with the finger, enabling easy positioning of the preparation.

The stage has a fixed condenser on its base.

The clamps may be replaced by a mechanical stage, part no.: 574200.

4.6 Focussing and stop plate

Focusing is done by two buttons controlled equally on one side or the other of the column, in quite a low position so that the hands remain resting on the table. It is the stage that moves vertically, without causing any inconvenience to the observer. Furthermore, its plane is horizontal: we can thus observe preparations in a liquid medium.

The pair of top buttons attached to the gearwheel driving the rack is used for coarse adjustment. It is fitted with a friction device whose release at the end of travel prevents any damage to the gears.

The travel of the stage upwards may be limited by an adjustable stop plate whose vertical control is visible between the back of the stage and the column.

To adjust it, the stop screw and locknut assembly are unscrewed. Place a standard preparation on the stage placed in the low position beforehand, then, by looking at the plane of the preparation from the side, raise the stage by moving it rapidly until the slide is almost in contact with the 40 x telescopic objective. Then tighten the stop screw fully and lock it with the lock nut. It then becomes impossible to move the stage further up with the coarse adjustment.

The pair of lower buttons operates the fine adjustment designed to perfect the focus. The objectives being parfocal, the buttons are sufficient to restore sharpness when the magnification is changed.

4.7 Condenser

The microscope is equipped under the stage with a fixed condenser, with an iris aperture and filter-holder. Its F-number: 0.65 corresponds to that of the 40x objective. It is therefore advisable to close the iris when using objectives with lower magnifications: 10x and 4x. Remember that:

- Closing the iris increases the contrast and the depth of field, but reduces the resolution
- The opening must not be greater than that of the objective
- For a preparation that is too light, the removable blue filter may be placed on the filter-holder.

5 Maintenance

5.1 Protection against dust

Immediately after use, cover the device with its dust cover. If dust has accumulated on the lenses, blow with a syringe or a rubber bulb, or use a sable brush: the lenses are fragile and can be scratched by dust.

5.2 Finger marks

Do not touch the lenses. If, however, finger marks are seen on them, remove them with a soft cloth dipped in xylol.

5.3 Humidity

Humidity, combined with temperature, is particularly harmful for mounted optics. It is advisable to store the microscopes in a dry area, and in a particularly open climate, in their polystyrene cabinets along with a moisture-absorptive pouch.

5.4 Dismantling


A microscope is a precision instrument that is delicately manufactured. If it needs to be dismantled (maintenance), the instrument must be returned to the supplier or to an appropriate workshop.

6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*