



Fixed Speed Vortex Mixer

MX-E



Specifications

Items	Specifications
Voltage [VAC]	100-240
Frequency [Hz]	50/60
Power [W]	20
Shaking movement	Orbital
Orbital diameter [mm]	4.8
Motor rating input [W]	12
Motor rating output [W]	8
Speed range [rpm]	3000
Run type	Touch operation
Dimensions (mm)	134x134x70.5
Weight [kg]	0.8
Permitted ambient temperature[°C]	5-40
Permitted relative humidity	80%
Protection class acc. to DIN 60529	IP21

Please read carefully before use, and follow all operating and safety instructions!

Operation

- Remove the transportation lock, when the first time use.
- Put the instrument on a stable and plug in the main power.

Transportation lock



- If a test tube is pressed into the mixing head vertically, the instrument will begin to work.
- If the pressure force removed, the instrument will stop to work.
- Disconnect the power.

Safety Precautions

The detailed precautions described in the User Manual are carefully developed in an attempt to cover all the possible risks. Please read carefully before operating this centrifuge.

- Ensure that only trained staff works with the instrument.
- Never use the centrifuge in any manner not specified in these instructions.
- If any abnormal condition occurs during operation, please stop it immediately and contact our service center.

Warranty

This centrifuge and rotors are guaranteed for one year from the date of delivery provided that they have been operated and maintained properly.

For all repairs, adjustments, or replacement parts, please contact:

SAV – Le Laborantin
468 rue Jacques Monod
BP 1900
27019 EVREUX CEDEX France

0 8256 888 111*

** 0,15 € TTC/min. + prix appel*

12300423

Distribué par
Le Laborantin
www.lelaborantin.com
clients@lelaborantin.com